

Customer Report

Wednesday, October 07, 2015

C
O
N
T
A
C
T

Viroblock

Chemin des Aulx 18.

Plans-les-Ouates

CH-1228

Switzerland

Thierry Pelet

41 22 884 83 44

t.pelet@viroblock.com

Project Title

Antimicrobial Testing

ID

0915-BOY-01 -- 1

Entry Date

9/24/2015

Project Summary

The **AATCC TM 100** test method is designed to measure the antimicrobial properties of textile or absorbent material incubated with selected microorganisms. The basis of the test methods is the incubation of the microorganism inoculum in contact with the test sample for a duration of up to 24 hours without drying. Following this exposure, the inoculated microorganisms are recovered and the concentration of the organisms is determined. *Candida albicans* was tested according to the standard method, culturing *C. albicans* prior to testing was conducted as required by the organism. The antimicrobial performance is determined by comparison of the recovered organisms from the test samples at time 0, and treated material after selected time points and is reported as a percent value relative to the control sample material.

Recommended Reading

Online Resource for Product Development, Testing, and Inquiry; *The Wily Microbe*

Guidance on antimicrobial preservation

<http://wily-microbe.situbiosciences.com/15-microbial-control/>
<http://wily-microbe.situbiosciences.com/34-microbial-control-testing/>

Antimicrobial testing with textiles

<http://wily-microbe.situbiosciences.com/280-textile-testing-antimicrobials/>

Sample List

Method Name

<i>Sample #</i>	<i>Sample Name</i>	<i>Sample Notes</i>
AATCC 100 Assessment of Antibacterial Finishes on Textile Materials		
1	TBD Fabric	Treated fabric in Triplicate
2	UTC	Control in Triplicate

Project - ATCC Ordering & Culturing

3 C. albicans adaptation

Result Table

Contact	Viroblock	Thierry Pelet	41 22 884 83 44
Title	Antimicrobial Testing		
Project ID	0915-BOY-01 -- 1	Entry Date 9/24/2015	Test Start Date 9/24/2015

Result Table *

Test Method AATCC 100 Assessment of Antibacterial Finishes on Textile Materials

Sample #	1	TBD Fabric	Interval	Result
Inoculum	C. albicans (10231)			
<i>Notes Section</i>				
replicate 1			24 hr	99.68 % Reduction
replicate 2			24 hr	99.91 % Reduction
replicate 3			24 hr	99.86 % Reduction
Inoculum	S. aureus (6538)			
<i>Notes Section</i>				
replicate 1			24 hr	99.99 % Reduction
replicate 2			24 hr	99.99 % Reduction
replicate 3			24 hr	99.99 % Reduction

Sample #	2	UTC	Interval	Result
Inoculum	C. albicans (10231)			
<i>Notes Section</i>				
replicate 1			0 hr	520000 CFU/ml
replicate 2			0 hr	480000 CFU/ml
replicate 3			0 hr	460000 CFU/ml
replicate 1			24 hr	920000 CFU/ml
replicate 2			24 hr	970000 CFU/ml
replicate 3			24 hr	1100000 CFU/ml
Inoculum	S. aureus (6538)			
<i>Notes Section</i>				
replicate 1			0 hr	480000 CFU/ml
replicate 2			0 hr	490000 CFU/ml
replicate 3			0 hr	510000 CFU/ml
replicate 1			24 hr	11000000 CFU/ml
replicate 2			24 hr	12000000 CFU/ml
replicate 3			24 hr	11000000 CFU/ml

Result Table *

Test Method Project - ATCC Ordering & Culturing

Sample # 3 C. albicans adaptation

Interval

Result

Inoculum None ()

Notes Section

complete

Result Table *

Test Method - Additional Information

AATCC 100 Assessment of Antibacterial Finishes on Textile Materials

Method Conventions

Test results are reported as the % Reduction of the test sample microorganism counts relative to the control sample microorganism counts. Samples demonstrating no reduction are reported as "0% Reduction". Microorganism concentrations are reported as the number of colony forming units (cfu) per swatch of material. The number of material swatches used in each test is reported. Standard neutralizing media used is 1x or 1:10x D/E Neutralization Broth. Modification to the test media may be done to facilitate the use of alternative microorganisms in the test protocol.

Terminology

activity, n..of an antibacterial agent, a measure of effectiveness of the agent.

antibacterial agent, n..in textiles, any chemical which kills bacteria (bactericide) or interferes with microorganism growth.

Percent Reduction (% Reduction) – the measure of the difference between the microorganisms counted for control and test samples during a specific test period. Testing is conducted using identical conditions unless noted for the samples types. Concentrations are determined as CFU/ml.

Colony forming units (cfu)/ volume (milliliter – ml) : standard units of microorganisms concentration.

CFU/swatch - colony forming units / swatch of material tested

Untreated Control (UTC) - untreated control sample material used to demonstrate normal test performance, showing robust microorganism growth.

Interval - represents the point or time point from which the result value was determine; T0 indicates that the result is from the soonest possible time from inoculation to recovery of the inoculated sample (typically < 5min).

Result - the result is the measure of change or abundance. Result units indicate the actual measurements, frequently relative to a control value depending on the method or test requirements.

Tests of uncertainty for microbiological methods are based on empirically determined data either from replicate tests conducted in the laboratory or provided by external sources such as Interlaboratory Studies. When possible, external references will be used to guide the test method uncertainty estimate.

Test measurement uncertainty is based on the established standard provided in SOP 050410 Measurement Uncertainty Estimates.

Measurement Uncertainty was determined by the method of reproducibility of replicates (using T0 and T1 and/or T24 hour samples) (See A2LA Reference G108)

Expanded Uncertainty for the test method of k=2 is for a 95% confidence of Log10 (0.034).

Uncertainty Values in CFU are obtained by converting CFU counts (C1) to Log10 values (Log10 C1); multiply this result by the Expanded Uncertainty (Log10 C1*EU), then add and subtract from the Log10 value (Log C1 +/- (Log10 C1*EU)); convert back by taking the anti-log ($1 \times 10^{\text{Log C1} \pm (\text{Log10 C1} \times \text{EU})}$) which provides the upper and lower limits of 95% confidence in CFU.

AATCC provides a statement of precision and bias as follows:

Studies (see 13.9) indicate the following within-laboratory precision of the Standard Plate Count (SPC) Test: (a) among-analyst variation of 18% and (b) within-analyst variation of 8%. REF (Peeler, J. T.; Leslie, J. W.; Messer, J. W. *Replicate counting errors by analysts and bacterial colony counters. J. Food Protection, Vol. 45, 1982, pp 238-240.*)

* This report is governed by and incorporates by reference, the conditions of testing as posted on the date of issuance and is intended for your exclusive use. Any Copying or replication of this report to or for any other person or entity, or use of our company name or Service Mark is permitted only with our prior written consent. This report sets forth our findings solely with respect to test samples identified herein. The results set forth in this report are not indicative or representative of the quality or characteristics of the lot from which a test sample was taken or any similar identical product unless specifically and expressly noted. Our report includes all tests requested and the results there of based upon the information provided. You have 60 days from the date of issuance of this report to notify us of any material error or omission caused by our handling of the samples, provided however that such notice shall be in writing and shall specifically address the issue you wish to raise. A failure to raise such issue within the prescribed time shall constitute your unqualified acceptance of the completeness of this report, the test conducted and the correctness of the report contents.

d.p satchell Ph.D.
Technology Director